Technical Data Sheet



NOROX®MCP-75

Methyl ethyl ketone peroxide; Cumylhydroperoxide CAS#1338-23-4; 80-15-9 Liquid mixture

Description

NOROX®MCP-75 is a solution of Methyl ethyl ketone peroxides and Cumyl hydroperoxide in phlegmatizer. NOROX®MCP-75 is a liquid polymerization initiator for the room temperature cure of unsaturated polyester and vinyl ester resins with the following advantages/properties:

- Very low peak exotherm temperature
- Slower gel time than standard MEKP
- Excellent final cure both in thin and thick laminates
- Low impurity levels (water, MEK, salts)
- Less shrinkage and stress problems

Technical Data

Appearance	Pale yellow liquid		
Active oxygen	8.7 – 9.0%		
Density at 20°C	1.09 – 1.13 g/cm ³		
Flash point (Seta closed cup)	> 65°C (149°F)		
Critical temperature (SADT)	60°C (140°F)		
Soluble in	Oxygenated organic solvents		
Slightly soluble in	Water		
Recommended storage temperature	Max 30°C.		
Maintenance of activity at 30°C as from date of delivery	6 months		

Application

NOROX®MCP-75 has a flat and low exotherm curve which ends up in a good 24 hour cure. This curing behavior reduces stress, cracking and shrinkage problems in laminates & castings. Less shrinkage also means less fibre print through of gelcoated surfaces. The gel time and cure time is slower compared to standard MEKP products. NOROX®MCP-75 can be used when you want to build up thick laminates in one step, or if the laminate has big variations in thickness, or when you have high resin loadings. Typical applications are vacuum bag moulding of big units, filament winding of pipes and tanks with thicker laminates or polyester concrete systems with high resin-to-filler ratio.

Technical Data Sheet




Measurements

A reactivity test with an unsaturated polyester resin gave the following results:

Resin:	Orthophthalic polyester	Temperature:		21°C	
Initiator %:1.0		Accelerator %		1% (1% cobalt)	
Initiat	tor	Gel	Time to		Peak
		time	peak		exotherm temp
		min	Min		°C
NOR	OX®MCP-75	22	69		31
NOR	OX [®] KP-9	11	48	•	47
NOR	OX [®] MEKP-925	15	59		43

Reactivity of NOROX® Initiators:



Standard Packaging

The standard package sizes of NOROX®MCP-75 are 5 kg (8 lb) and 25 kg (40 lb) polyethylene bottles.

Disclaimer

This information and all further technical advice are reflecting our present knowledge and experience based on internal tests with local raw materials with the purpose to inform about our products and applications. The information should not be construed as guaranteeing specific properties of products described or their suitability for a particular application, nor as providing complete instructions for use. The information implies no guarantee for product and shelf life properties, nor any liability or other legal responsibility on our part, including with regard to existing third party intellectual property rights, especially patent rights. We reserve the right to make any changes according to technological progress or further developments.

Application and usage of our products based on our technical advice is out of our control and sole responsibility of the user. The user is not released from the obligation to conduct careful inspection and testing of incoming goods in order to verify the suitability for the intended application.

United Initiators **EU**T: +49 89 74422 237
F: +49 89 74422 6237
cs-initiators.eu@united-in.com

United Initiators
Nafta
T: +1 800 231 2702
F: +1 440 323 0898
cs-initiators.nafta@united-in.com

United Initiators
China
T: +86 20 6131 1370
F: +86 139 2503 8952
cs-initiators.cn@united-in.com

United Initiators **Australia**T: +61 2 9316 0046
F: +61 2 9316 0034
cs-initiators.au@united-in.com

www.united-initiators.com

Revision number: 1.0. Date: 03.03.14. Device M: TDS.